

Product Introduction

Model : CCE-3306B

Product Description



Product Features :

1. Portable and delicate design, easy to operate
2. Accurate pulse pace identification function
3. High accuracy digital filter, automatic baseline adjustment
4. Working modes: Manual, Automatic, Arrhythmia Analysis, Storage
5. 1120mm, six channel format recording, excellent automatic interpretation
6. 800x480 graphic 7 inch color LCD to simultaneously display ECG information
7. 40 patient cases storage and replay
8. Detailed patient information record
9. Adapt to 110-230V,50/60Hz power supply. Built-in rechargeable Li-ion battery
10. RS232 and USB interface, support for PC ECG software function. (optional)

Technical Specification :

| | | | |
|-------------------------|--|----------------------------|---|
| Lead | Standard 12 leads | Time constant | $\geq 3.2s$ |
| Lead acquisition | 12bit/1000Hz (synchronously 12 leads) | Frequency response | 0.05~160Hz(-3db) |
| Work mode | Man / Auto/Analysis/Storage | Noise level | $< 15\mu Vp-p$ |
| Filter | AC Filter : 50Hz / 60Hz EMG Filter : 25Hz / 45Hz Anti-Drift Filter : 0.15Hz (Adaptive) | Inter-channel Interference | $\leq 0.5mm$ |
| CMMR | $\geq 120dB$ (with AC FILTER) | Sensitivity | 5, 10, 20, 40, Auto (mm/mV) |
| Input circuit | Floating; Protection circuit against Defibrillator effect | Recording mode | 3ch, 3ch+, 3ch+++, 6ch Default Man. Mode is 3ch format Default Auto. Mode is 6ch format |
| Input Impedance | $\geq 50M\Omega$ | Paper speed | 6.25、12.5、25、50mm/s ($\pm 3\%$) |
| Input circuit current | $\leq 0.05\mu A$ | Paper | 112mm *20m high resolution thermal roll paper |
| Patient current leakage | $< 10\mu A$ | Communication | RS232/ USB port |
| Calibrating voltage | 1mV $\pm 2\%$ | Power Supply | AC 110-230V($\pm 10\%$),50/60Hz ($\pm 1Hz$),60VA DC Built-in rechargeable Li-ion battery,14.4V (2000mAh) |
| Voltage tolerance | $\geq \pm 500mV$ | Dimension/Weight | 300mm \times 230mm \times 75mm / 2.8Kg |

Packing List :

| |
|--|
| 1. Main Unit : 1 set |
| 2. Patient Cable : 1pc |
| 3. Chest Electrodes : 1set (6 pcs) |
| 4. Limb Electrodes : 1set (4 pcs) |
| 5. Power Cable : 1pc |
| 6. Grounding Cable : 1pc |
| 7. Fuse : 2pcs |
| 8. Paper Axis : 3pcs |
| 9. Thermal Paper : 1 roll |
| 10. Use Manual : 1 copy |
| 11. The factory Inspection Report : 1 copy |

Picture for Accessories :

